



DeltaT comparison

PC software for your Controlled Rate Freezer

- Unlimited freezing profiles
- Graphical data in real time
- Graphical data after the run
- Traceability of records
- Password protection
- Pre and post run annotation



DeltaT provides a PC interface for the Planer range of controlled rate freezers allowing you to rapidly develop and maintain profiles on your PC and view freezing progress in real time.

New to Planer is DeltaT v7, designed today's needs, with batch mode and networking capabilities.

DeltaT can be used to store files (only limited by the size of the hard drive) on your PC and to print using a standard printer for a manual record.

The full version of DeltaT v6 is provided with the Kryo 550, 750 and 1060 machines.

DeltaT-lite is available free of charge from www.planer.com.

DeltaT v6 is only available via download from www.planer.com. Once downloaded, you need to email licence@planer.com with the full details from the DeltaT licensing menu, to have the software licensed.

DeltaT v7 demo version is available via download from www.planer.com. Once downloaded, you need to email licence@planer.com with the full details from the DeltaT licensing menu, to have the software licensed.

DeltaT v7 can be used to view a live graph via a network (LAN) designed just to send data from the CRF by email.

DeltaT v7 can let the user know if, for example, the freezer has gone into alarm mode, and also how a run is progressing, emailing the run data after a completed run.

DeltaT is available in four versions shown in the table below

	DeltaT- lite v6	DeltaT v6	DeltaTiQ v6 (Distributors only)	DeltaT v7
Create profiles	yes	yes	yes	yes
Run Profiles	yes	yes	yes	yes
View kryofiles	yes	yes	yes	yes
Print kryofiles	x	yes	yes	yes
Export kryofiles	x	yes	yes	yes
Calibration	x	x	yes	yes
Run qualification	x	x	yes	yes
Batch mode	x	x	x	yes
Network capabilities	x	x	x	yes
Compatibility with the MR6 controller	yes	yes	yes	yes