


# MVE High Efficiency 800 -190°C Series



## Features include


- Dry sample storage
- -190°C top box temperature
- Lowest liftover height
- Largest LN<sub>2</sub> capacity at vapor platform
- Optional Battery Backup



- Optional One-Tier and Two-Tier Steps

## Featuring the TEC 3000

The MVE High Efficiency 800 Freezer Series provides vapor storage temperature at -190°C. Available in three unique sizes these freezers provide maximum storage density and provide the industry's longest hold time.



## Temperature Test Graph

\* Temp Test indicates typical performance of 800 Series -190°C freezer with full inventory system and factory recommended level settings. Actual performance may vary with atmospheric conditions and usage.

[www.chartbiomed.com](http://www.chartbiomed.com)  
844-MVE-CRYO



	MVE 815P-190			MVE 818P-190			MVE 819P-190		
Maximum Storage Capacities									
1.2 & 2 ml Vials (Internally Threaded)	15,600			18,200			19,500		
Quantity of Large Racks 100 cell boxes	12			12			12		
Quantity of Mini Racks 25 cell boxes	4			4			4		
Number of Shelves per Rack	12			14			15		
Performance									
LN <sub>2</sub> Capacity w/o Inventory L est.	370			420			463		
LN <sub>2</sub> Capacity at Vapor Platform L est.	52			55			55		
Unit Dimensions									
Neck Opening in. (mm)	12.5 (317)			12.5 (317)			12.5 (317)		
Usable Internal Height in. (mm)	26.5 (673)			30.7 (781)			34.5 (877)		
Inner Diameter in. (mm)	28.8 (731)			28.8 (731)			28.8 (731)		
Overall Height w/Autofill in. (mm)	57.0 (1448)			61.3 (1556)			65.0 (1651)		
Door Width Requirement** in. (mm)	32.0 (813)			32.0 (813)			32.0 (813)		
Door Width Requirement** with handles in. (mm)	33.2 (843)			33.2 (843)			33.2 (843)		
Weight Empty* lb. (kg)	475 (215)			495 (225)			515 (234)		
Weight Liquid Full* lb. (kg) est.	1134 (514)			1168 (530)			1340 (608)		
Blood Bag Capacities	Total Bags	Bags/ Frame	No. Frames	Total Bags	Bags/ Frame	No. Frames	Total Bags	Bags/ Frame	No. Frames
791 OS/U Medsep (25 ml)	1,224	6	204	1,632	8	204	1,836	9	204
4R9951 (50 ml)	768	6	128	896	7	128	1024	8	128
4R9953 (250 ml)	416	4	104	416	4	104	520	5	104
4R9955 (500 ml)	304	4	76	304	4	76	380	5	76
DF200 (200 ml)	236	4	59	236	4	59	295	5	59

**TWO Year Parts Warranty • FIVE Year Vacuum Warranty**


Conforms to MDD 93/42/EEC, the Medical Device Directive for the EU. Freezer systems UL/C-UL Listed.

\* Without inventory



\*\*Minimum width required for vessel to pass through opening. Footprint may vary. Contact Tech Service for detailed drawings.

## Rack Layout

Square Rack Layout (P)



4R9951 Layout (R)



Visit [www.chartbiomed.com](http://www.chartbiomed.com) for more information.

© 2019 Chart Inc. reserves the right to discontinue its products, or change the prices, materials, equipment, quality, descriptions, specifications and/or processes to its products at any time without prior notice and with no further obligation or consequence. All rights not expressly stated herein are reserved by us, as applicable.