

CERTIFICATE OF TREATMENT BY GAMMA RADIATION

WE UNDERSIGNED :

Synergy Health S.A.S.

MIN 712 - Les ARNAVAUX

13323 MARSEILLE CEDEX 14 - FRANCE



CERTIFY THAT WE TREATED BY GAMMA RADIATION ACCORDING TO THE SPECIFIC CUSTOMER'S REQUIREMENTS AND TO :

- the specification of treatment # 496604P
- the specifications of ISO 13485 and ISO 11137-1 standards
- the results of the dosemapping # 169936 of 07/09/2011

THE FOLLOWING PRODUCTS : (according to the customer's indication)

Customer	: PROMEPLA
Product	: DEVICE OOG00199SR148
Customer's reference	: ORDER # 99/PLANER OF 04/04/2013
Quantity	: 3 PALLETS
Irradiation date	: 2013.04.06
Irradiation dose	: 27.7 kGy to 37.6 kGy
Irradiation batch number	: 19028201

The control of the applied radiation dose is done by Synergy Health using Red Perspex dosimeters calibrated by the English National Physical Laboratory.

Synergy Health Marseille,

A. NOEL

Process Control Officer

A handwritten signature in blue ink, appearing to read 'A. Noel', is written over the printed name and title.

Certificate # 19028201 / 1

S. LE GONIDEC

Quality Assistant

A handwritten signature in blue ink, appearing to read 'S. Le Gonidec', is written over the printed name and title.