



<http://www.steris-ast.com>

Certificate of Irradiation

Date Issued: 22-Feb-2021

FR01S12538885-1-1

This is to certify that Synergy Health Marseille, a STERIS Company has where appropriate delivered an irradiation process in accordance with the current certified standards:

EN ISO 11137-1 Sterilisation of Health Care Products
EN ISO 13485 Quality System - Medical Devices

PROMEPLA
Le Copori
9 avenue ALbert II
MONACO
98000
MONACO

Order Information

Account Number:	104501
Synergy Health Sales Part Reference:	1129014
Customer Reference Number:	013A-21MG - 13/01/2021
Product Description:	ESSAI CARTO - PROMEPLA
Quantity Received:	1
Customer Minimum Specification kGy:	25.0
Customer Maximum Specification kGy:	45.0
Customer Unit Lot/Batch Number:	G005194

Irradiation Data

Date and Time of Irradiation:	22-FEB-2021 02:15
Calculated Minimum Dose kGy:	28.2
Calculated Maximum Dose kGy:	39.8

CN200115
BOTTLE ASSY WITH LABEL (BT37 MK2)
Supp Ref:OOG00799SR148
Bin:QS/1 GRN:GR1078

Qty:1,392 Each

Purchase Order
#PN188
MTMS PO #
PN045647

Irradiation Release Authorised By Synergy Health Marseille SAS, a STERIS Company

Processing Site: M.I.N. 712 - Arnavaux, , Marseille Cedex 14, 13323 Phone No: +33 (0) 4 91 214 214

Registered Office: M.I.N. 712 - Arnavaux, 13323 Marseille Cedex 14, FRANCE
N° TVA: FR59 343 092 540

PLANER

PLANER Ltd, 110 Windmill Road, Sunbury-on-Thames
Middlesex, TW16 7HD, United Kingdom
Tel : +44(0)1932 755000 Fax : +44(0)1932 755001
E-mail : sales@planer.co.uk Web : www.planer.com



DESCRIPTION : STERILIZED BOTTLE ASSEMBLY- SINGLE TUBE VERSION

REF CN200115

RX ONLY

CATALOGUE NUMBER

LOT G005194
BATCH CODE

STERILE




STERILIZED USING IRRADIATION

DO NOT REUSE

 **2023-05-28**
USE BY



STORE IN A DARK PLACE AT ROOM TEMPERATURE



Do not use if package is damaged.
Nepoužívejte, je-li balení poškozeno.
Må ikke anvendes, hvis emballagen er beskadiget.
Niet gebruiken indien de verpakking beschadigd is.
Ei saa käyttää, jos pakkaus on vahingoittunut.
Ne pas utiliser si l'emballage est endommagé.
Inhalt bei beschädigter Verpackung nicht verwenden.



Μην το χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά.
Non utilizzare se la confezione è danneggiata.
Nie stosować, jeżeli opakowanie zostało uszkodzone.
Não utilizar se a embalagem estiver danificada.
No utilizar si el envase está dañado.
Får inte användas om förpackningen är skadad.
Nepoužívať, ak je balenie poškodené.