



<http://www.steris-ast.com>

Certificate of Irradiation

Date Issued: 28-Jul-2021

FR01S12591305-1-1

This is to certify that Synergy Health Marseille, a STERIS Company has where appropriate delivered an irradiation process in accordance with the current certified standards:

EN ISO 11137-1 Sterilisation of Health Care Products
EN ISO 13485 Quality System - Medical Devices

PROMEPLA
Le Copori
9 avenue ALbert II
MONACO
98000
MONACO

Order Information

Account Number:	104501
Synergy Health Sales Part Reference:	1134451
Customer Reference Number:	433/PLANER PLC - 22/07/2021
Product Description:	496633P - OOG00799SR148
Validation Reference:	S12538885
Quantity Received:	7
Customer Minimum Specification kGy:	25.0
Customer Maximum Specification kGy:	45.0
Customer Unit Lot/Batch Number:	G005724

Irradiation Data

Date and Time of Irradiation:	28-JUL-2021 09:48
Calculated Minimum Dose kGy:	26.2
Calculated Maximum Dose kGy:	43.9

Irradiation Release Authorised By Synergy Health Marseille SAS, a STERIS Company

Processing Site: M.I.N. 712 - Arnavaux, , Marseille Cedex 14, 13323 Phone No: +33 (0) 4 91 214 214

Registered Office: M.I.N. 712 - Arnavaux, 13323 Marseille Cedex 14, FRANCE

N° TVA: FR59 343 092 540

Détail de la commande: 433

Référence	Ref client	Quantité	Nb carton	Num lot	N° palette
CLIENT: PLANER PLC					
OOG00799SR148		480	10	G005724	R210091
OOG00799SR148		480	10	G005724	R210092
OOG00799SR148		480	10	G005724	R210093
OOG00799SR148		480	10	G005724	R210094
OOG00799SR148		480	10	G005724	R210095
OOG00799SR148		480	10	G005724	R210096
OOG00799SR148		480	10	G005724	R210097
		3360	70		

PROMEPLA S.A.M.

Monaco le 18/08/21

9, Avenue Albert II
98000 Monaco

TEL : 00 377 97 98 42 32
FAX : 00 377 92 05 61 50
E-mail : sales@promempla.com
Site : www.promempla.com

Packing List N° BL211859

PLANER LTD
110 WINDMILL ROAD
MIDDLESEX
TW16 7HD SUNBURY-ON-THAMES
ANGLETERRE (UK)

Delivery address

PLANER PLC
110 WINDMILL ROAD
UNIT F, DOLPHIN IND ESTATE
TW16 7HD SUNBURY-ON-THAMES
ANGLETERRE (UK)

Your Ref : PN046271 MAY Code : CPLANER

Shipment : MONACO LOGISTIQUE
N° Tracking : 80cmX1m20X1m94
Number of Cases : 70
Handling Unit : 7
Total Weight : 700

Part N°	Y/ Part N°	Description	Batch N° Compl.	Batch N°	Qty	Still to deliver
0000799SR148	CN200115	BT37 MKII - BOTTLE ASSEMBLY CUSTOMS 901839000/ EORI NUMBER: FR00 00 0905900010/ FRANCE ORIGIN		G005724	3360	

PROMEPLA S.A.M.

9, Avenue Albert II
Immeuble Le Copert

98000 MONACO

Phone number : (+377) 97 98 42 32

Mail : contact@promempla.com

CERTIFICATE OF COMPLIANCE

Bearing all circumstances beyond our control, we hereby certify that all the above-mentioned products bearing a batch number are, in conformity with the technical specifications and have not undergone any modifications whatsoever with respect to raw materials, processing or technical characteristics and are in accordance with our NF EN ISO13485 : 2016 Quality System.

EXTRACT OF GENERAL TERMS OF SALES

All our goods travel at the address's own risks.

S.A. Montégasque au capital de 588.420 € - N.I.E : FR290000009059 - S.S.E.E. 4669 B0 1526 - R.C.I 69 S 1239 -

PLANER

PLANER Ltd, 110 Windmill Road, Sunbury-on-Thames
Middlesex, TW16 7HD, United Kingdom

Tel : +44(0)1932 755000

E-mail : sales@planer.com

Web : www.planer.com



DESCRIPTION : STERILIZED BOTTLE ASSEMBLY - SINGLE TUBE VERSION

REF CN200115

RX ONLY



1639

LOT G005724

BATCH CODE

2021-05-28



DATE OF

MANUFACTURE

2023-10-28



USE BY

STERILE R

STERILIZED USING IRRADIATION DO NOT REUSE



DO NOT USE IF THE PACKAGE IS DAMAGED



STORE IN A DARK PLACE AT ROOM TEMPERATURE

Do not use if package is damaged.

Nepoužívejte, je-li balení poškozeno.

Mã ikke anvendes, hvis emballagen er beskadiget.

Niet gebruiken indien de verpakking beschadigd is.

Ei saa käyttää jos pakkaus on vaurioitunut.

Ne pas utiliser si l'emballage est endommagé.

Inhalt bei beschädigter Verpackung nicht verwenden.

Μην το χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιώ.

Non utilizzare se la confezione è danneggiata.

Nie stosować, jeżeli opakowanie zostało uszkodzone.

Não utilizar se a embalagem estiver danificada.

No utilizar si el envase está dañado.

Får inte användas om förpackningen är skadad.

Nepoužívat, ak je balenie poškodené.

EC REP

Advent Ltd, Lower Business centre, 2nd Flr.,
Tower Street, Swatar, BKR 4013 Malta

GRN: GR2182

Supplier	PA169 Promepla Sa	Receipt date	27/08/2021	Created From	Purchase #PN238	Order
Supplier advice note	BL211859	Department				

PO line#	Count items ?	Location	Part#	Store	Bin	Serial/Lot#	Expiration	Qty
1	No	Production	CN200115	QS	QS/1			3,360
		BOTTLE ASSY WITH LABEL (BT37 MK2)						